

TWIK3 Antibody

Catalog No: #47968



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TWIK3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	Hu,Ms,Rt
Specificity	Recognizes endogenous levels of TWIK3 protein.
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human TWIK3.
Target Name	KCNK7
Other Names	Potassium channel subfamily K member 7
Accession No.	Swiss-Prot#:Q9Y2U2NCBI Gene ID:10089
Uniprot	Q9Y2U2
GeneID	10089;
Calculated MW	36KD
Concentration	1 mg/ml
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of TWIK3 expression in HCT116 (A), A549 (B), H9C2 (C), AML12 (D) whole cell lysates

Note: This product is for in vitro research use only