

Ribosomal Protein S9 Polyclonal Antibody

Catalog No: #41412

Package Size: #41412-1 50ul #41412-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribosomal Protein S9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Hu Ms Rt
Specificity	Ribosomal Protein S9 Polyclonal Antibody detects endogenous levels of Ribosomal Protein S9 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ribosomal Protein S9.
Target Name	Ribosomal Protein S9
Other Names	RPS9; 40S ribosomal protein S9
Accession No.	Swiss-Prot: P46781NCBI Gene ID: 6203
Uniprot	P46781
GeneID	6203;
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

Images



Western Blot analysis of HuvEc cells using Ribosomal Protein S9 Polyclonal Antibody

Note: This product is for in vitro research use only