

Human IgG heavy chain Antibody

Catalog No: #48056

Package Size: #48056-1 50ul #48056-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Human IgG heavy chain Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purified
Applications	WB
Species Reactivity	Hu
Immunogen Description	Human IgG
Accession No.	Swiss-Prot#:
Calculated MW	50 kDa
Formulation	1*TBS (pH7.4), 0.5%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:5,000

Images



Western blot analysis on human IgG protein using anti-Human IgG heavy chain Mouse mAb (Cat. # M0806-1).

References

Note: This product is for in vitro research use only