

RyR-2 Polyclonal Antibody

Catalog No: #41420

Package Size: #41420-1 50ul #41420-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RyR-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	RyR-2 Polyclonal Antibody detects endogenous levels of RyR-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human RyR-2 around the non-phosphorylation site of S2808.
Target Name	RyR-2
Other Names	RYR2; Ryanodine receptor 2; RYR-2; RyR2; hRYR-2; Cardiac muscle ryanodine receptor; Cardiac muscle ryanodine receptor-calcium release channel; Type 2 ryanodine receptor
Accession No.	Swiss-Prot: Q92736NCBI Gene ID: 6262
Uniprot	Q92736
GeneID	6262;
SDS-PAGE MW	200-300kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1:500-2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

Images



Western Blot analysis of HepG2-UV MOUSE-BRAIN AD293
PC-3 cells using RyR-2 Polyclonal Antibody

Note: This product is for in vitro research use only