

Her2/ErbB2 Antibody

Catalog No: #48245



Package Size: #48245-1 50ul #48245-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Her2/ErbB2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Peptide affinity purified
Applications	WB, ICC, IHC, FC
Species Reactivity	Hu
Immunogen Description	peptide
Other Names	Verb b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog antibody C erb B2/neu protein antibody CD340 antibody CD340 antigen antibody Cerb B2/neu protein antibody CerbB2 antibody Erb b2 receptor tyrosine kinase 2 antibody ERBB2 antibody ERBB2_HUMAN antibody HER 2 antibody HER 2/NEU antibody HER2 antibody Herstatin antibody Human epidermal growth factor receptor 2 antibody Metastatic lymph node gene 19 protein antibody MLN 19 antibody MLN19 antibody NEU antibody NEU proto oncogene antibody Neuro/glioblastoma derived oncogene homolog antibody Neuroblastoma/glioblastoma derived oncogene homolog antibody NGL antibody p185erbB2 antibody Proto-oncogene c-ErbB-2 antibody Proto-oncogene Neu antibody Receptor tyrosine-protein kinase erbB-2 antibody TKR1 antibody Tyrosine kinase type cell surface receptor HER2 antibody Tyrosine kinase-type cell surface receptor HER2 antibody V erb b2 avian erythroblastic leukemia viral oncogene homolog 2 (neuro/glioblastoma derived oncogene homolog) antibody V erb b2 avian erythroblastic leukemia viral oncogene homolog 2 antibody V erb b2 avian erythroblastic leukemia viral oncoprotein 2 antibody V erb b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian) antibody V erb b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog antibody Verb b2 erythroblastic leukemia viral oncogene homolog 2, neuro/glioblastoma derived oncogene homolog (avian) antibody
Accession No.	Swiss-Prot#:P04626
Uniprot	P04626
GeneID	2064;
Formulation	1*TBS (pH7.4), 1%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

References

Note: This product is for in vitro research use only