

Synaptotagmin Polyclonal Antibody

Catalog No: #41471



Package Size: #41471-1 50ul #41471-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Synaptotagmin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Synaptotagmin Polyclonal Antibody detects endogenous levels of Synaptotagmin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Synaptotagmin 1/2 around the non-phosphorylation site of S309/306.
Target Name	Synaptotagmin
Other Names	SYT1; SVP65; SYT; Synaptotagmin-1; Synaptotagmin I; SytI; p65; SYT2; Synaptotagmin-2; Synaptotagmin II; SytII
Accession No.	Swiss-Prot: P21579/Q8N9I0NCBI Gene ID: 6857/127833
Uniprot	P21579
GeneID	6857;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Synaptotagmin Polyclonal Antibody

Note: This product is for in vitro research use only