

TGFβ RII Polyclonal Antibody

Catalog No: #41492

Package Size: #41492-1 50ul #41492-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TGFβ RII Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	TGFβ RII Polyclonal Antibody detects endogenous levels of TGFβ RII protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human TGFβ RII around the non-phosphorylation site of S225.
Target Name	TGFβ RII
Other Names	TGFBR2; TGF-beta receptor type-2; TGFR-2; TGF-beta type II receptor; Transforming growth factor-beta receptor type II; TGF-beta receptor type II; TbetaR-II
Accession No.	Swiss-Prot: P37173NCBI Gene ID: 7048
Uniprot	P37173
GeneID	7048;
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of NIH-3T3 cells using TGF β RII Polyclonal Antibody

Note: This product is for in vitro research use only