

TNF- $\beta$  Polyclonal Antibody

Catalog No: #41505

Package Size: #41505-1 50ul #41505-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TNF- $\beta$ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	TNF- $\beta$ Polyclonal Antibody detects endogenous levels of TNF- $\beta$ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TNF- $\beta$ .
Target Name	TNF- $\beta$
Other Names	LTA; TNFB; TNFSF1; Lymphotoxin-alpha; LT-alpha; TNF-beta; Tumor necrosis factor ligand superfamily member 1
Accession No.	Swiss-Prot: P01374NCBI Gene ID: 4049
Uniprot	P01374
GeneID	4049;
SDS-PAGE MW	22kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

## Images

Western Blot analysis of COS7 cells using TNF- $\beta$  Polyclonal Antibody

---

Note: This product is for in vitro research use only