

APBB3 antibody

Catalog No: #22673

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	APBB3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 240 and 483 of Human APBB3
Target Name	APBB3
Accession No.	Swiss-Prot:O95704Gene ID:10307
Uniprot	O95704
GeneID	10307;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 52kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: A431
10% SDS PAGE
Primary antibody diluted at 1: 1000

Background

The protein encoded by this gene is a member of the APBB protein family. It is found in the cytoplasm and binds to the intracellular domain of the Alzheimer's disease beta-amyloid precursor protein (APP) as well as to other APP-like proteins. It is thought that the protein encoded by this gene may

modulate the internalization of APP. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq]

Note: This product is for in vitro research use only