

TP53INP1 Polyclonal Antibody

Catalog No: #41507

Package Size: #41507-1 50ul #41507-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TP53INP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	TP53INP1 Polyclonal Antibody detects endogenous levels of TP53INP1 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human TP53INP1.
Target Name	TP53INP1
Other Names	TP53INP1; P53DINP1; SIP; Tumor protein p53-inducible nuclear protein 1; Stress-induced protein; p53-dependent damage-inducible nuclear protein 1; p53DINP1
Accession No.	Swiss-Prot: Q96A56NCBI Gene ID: 94241
Uniprot	Q96A56
GeneID	94241;
SDS-PAGE MW	27kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using TP53INP1 Polyclonal Antibody

Note: This product is for in vitro research use only