

TRAF6 Polyclonal Antibody

Catalog No: #41510

Package Size: #41510-1 50ul #41510-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TRAF6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	TRAF6 Polyclonal Antibody detects endogenous levels of TRAF6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TRAF6.
Target Name	TRAF6
Other Names	TRAF6; RNF85; TNF receptor-associated factor 6; E3 ubiquitin-protein ligase TRAF6; Interleukin-1 signal transducer; RING finger protein 85
Accession No.	Swiss-Prot: Q9Y4K3NCBI Gene ID: 7189
Uniprot	Q9Y4K3
GeneID	7189;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using TRAF6 Polyclonal Antibody



Western Blot analysis of HeLa cells using TRAF6 Polyclonal Antibody

Note: This product is for in vitro research use only