

TRAP220 Polyclonal Antibody

Catalog No: #41513

Package Size: #41513-1 50ul #41513-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TRAP220 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	TRAP220 Polyclonal Antibody detects endogenous levels of TRAP220 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TRAP220.
Target Name	TRAP220
Other Names	MED1; ARC205; CRSP1; CRSP200; DRIP205; DRIP230; PBP; PPARBP; PPARGBP; RB18A; TRAP220; TRIP2; Mediator of RNA polymerase II transcription subunit 1; Activator-recruited cofactor 205 kDa component; ARC205; Mediator complex subunit 1; Peroxiso
Accession No.	Swiss-Prot: Q15648NCBI Gene ID: 5469
Uniprot	Q15648
GeneID	5469;
SDS-PAGE MW	170kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using TRAP220 Polyclonal Antibody

Note: This product is for in vitro research use only