

USP38 Polyclonal Antibody

Catalog No: #41525

Package Size: #41525-1 50ul #41525-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	USP38 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	USP38 Polyclonal Antibody detects endogenous levels of USP38 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human USP38.
Target Name	USP38
Other Names	USP38; KIAA1891; Ubiquitin carboxyl-terminal hydrolase 38; Deubiquitinating enzyme 38; HP43.8KD; Ubiquitin thioesterase 38; Ubiquitin-specific-processing protease 38
Accession No.	Swiss-Prot: Q8NB14NCBI Gene ID: 84640
Uniprot	Q8NB14
GeneID	84640;
SDS-PAGE MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using USP38 Polyclonal Antibody

Note: This product is for in vitro research use only