

UTP14A Polyclonal Antibody

Catalog No: #41527

Package Size: #41527-1 50ul #41527-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	UTP14A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	UTP14A Polyclonal Antibody detects endogenous levels of UTP14A protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human UTP14A.
Target Name	UTP14A
Other Names	UTP14A; SDCCAG16; U3 small nucleolar RNA-associated protein 14 homolog A; Antigen NY-CO-16; Serologically defined colon cancer antigen 16
Accession No.	Swiss-Prot: Q9BVJ6NCBI Gene ID: 10813
Uniprot	Q9BVJ6
GeneID	10813;
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using UTP14A Polyclonal Antibody

Note: This product is for in vitro research use only