

Wnt-1 Polyclonal Antibody

Catalog No: #41535

Package Size: #41535-1 50ul #41535-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Wnt-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms
Specificity	Wnt-1 Polyclonal Antibody detects endogenous levels of Wnt-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Wnt-1.
Target Name	Wnt-1
Other Names	WNT1; INT1; Proto-oncogene Wnt-1; Proto-oncogene Int-1 homolog
Accession No.	Swiss-Prot: P04628NCBI Gene ID: 7471
Uniprot	P04628
GeneID	7471;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using Wnt-1 Polyclonal Antibody



Western Blot analysis of NIH-3T3 HT29 cells using Wnt-1 Polyclonal Antibody

Note: This product is for in vitro research use only