

c-FLIP Polyclonal Antibody

Catalog No: #41569

Package Size: #41569-1 50ul #41569-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	c-FLIP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	c-FLIP Polyclonal Antibody detects endogenous levels of c-FLIP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human c-FLIP.
Target Name	c-FLIP
Other Names	CFLAR; CASH; CASP8AP1; CLARP; MRIT; CASP8 and FADD-like apoptosis regulator; Caspase homolog; CASH; Caspase-eight-related protein; Casper; Caspase-like apoptosis regulatory protein;
Accession No.	Swiss-Prot: O15519NCBI Gene ID: 8837
Uniprot	O15519
GeneID	8837;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of SW480 SKOV3 HepG2 cells using c-FLIP Polyclonal Antibody

Note: This product is for in vitro research use only