

FGF-9 Polyclonal Antibody

Catalog No: #41588

Package Size: #41588-1 50ul #41588-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FGF-9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	FGF-9 Polyclonal Antibody detects endogenous levels of FGF-9 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FGF-9.
Target Name	FGF-9
Other Names	FGF9; Fibroblast growth factor 9; FGF-9; Glia-activating factor; GAF; Heparin-binding growth factor 9; HBGF-9
Accession No.	Swiss-Prot: P31371NCBI Gene ID: 2254
Uniprot	P31371
GeneID	2254;
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of SW480 cells using FGF-9 Polyclonal Antibody

Note: This product is for in vitro research use only