

## FGF-5 Polyclonal Antibody

Catalog No: #41599

Package Size: #41599-1 50ul #41599-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | FGF-5 Polyclonal Antibody                                                                                             |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB ELISA                                                                                                              |
| Species Reactivity    | Hu                                                                                                                    |
| Specificity           | FGF-5 Polyclonal Antibody detects endogenous levels of FGF-5 protein.                                                 |
| Immunogen Type        | peptide                                                                                                               |
| Immunogen Description | Synthesized peptide derived from the C-terminal region of human FGF-5.                                                |
| Target Name           | FGF-5                                                                                                                 |
| Other Names           | FGF5; Fibroblast growth factor 5; FGF-5; Heparin-binding growth factor 5; HBGF-5; Smag-82                             |
| Accession No.         | Swiss-Prot: P12034NCBI Gene ID: 2250                                                                                  |
| Uniprot               | P12034                                                                                                                |
| GeneID                | 2250;                                                                                                                 |
| SDS-PAGE MW           | 30kd                                                                                                                  |
| Concentration         | 1mg/ml                                                                                                                |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of mouse liver cells using FGF-5 Polyclonal Antibody

---

Note: This product is for in vitro research use only