

TF Polyclonal Antibody

Catalog No: #41603

Package Size: #41603-1 50ul #41603-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	TF Polyclonal Antibody detects endogenous levels of TF protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TF.
Target Name	TF
Other Names	F3; Tissue factor; TF; Coagulation factor III; Thromboplastin; CD142
Accession No.	Swiss-Prot: P13726NCBI Gene ID: 2152
Uniprot	P13726
GeneID	2152;
SDS-PAGE MW	33kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of BXPC-3 mouse eye mouse lung cells using TF Polyclonal Antibody

Note: This product is for in vitro research use only