

HoxD3 Polyclonal Antibody

Catalog No: #41604

Package Size: #41604-1 50ul #41604-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HoxD3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	HoxD3 Polyclonal Antibody detects endogenous levels of HoxD3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human HoxD3.
Target Name	HoxD3
Other Names	HOXD3; HOX1D; HOX4A; Homeobox protein Hox-D3; Homeobox protein Hox-4A
Accession No.	Swiss-Prot: P31249NCBI Gene ID: 3232
Uniprot	P31249
GeneID	3232;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using HoxD3 Polyclonal Antibody

Note: This product is for in vitro research use only