

Cdk4 Polyclonal Antibody

Catalog No: #41607

Package Size: #41607-1 50ul #41607-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cdk4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Cdk4 Polyclonal Antibody detects endogenous levels of Cdk4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Cdk4.
Target Name	Cdk4
Other Names	CDK4; Cyclin-dependent kinase 4; Cell division protein kinase 4; PSK-J3
Accession No.	Swiss-Prot: P11802NCBI Gene ID: 1019
Uniprot	P11802
GeneID	1019;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 K562 cells using Cdk4 Polyclonal Antibody

Note: This product is for in vitro research use only