

Fascin 1 Polyclonal Antibody

Catalog No: #41620

Package Size: #41620-1 50ul #41620-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fascin 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Fascin 1 Polyclonal Antibody detects endogenous levels of Fascin 1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Fascin 1.
Target Name	Fascin 1
Other Names	FSCN1; FAN1; HSN; SNL; Fascin; 55 kDa actin-bundling protein; Singed-like protein; p55
Accession No.	Swiss-Prot: Q16658NCBI Gene ID: 6624
Uniprot	Q16658
GeneID	6624;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 HuvEc HT29 K562 L929 cells using Fascin 1 Polyclonal Antibody

Note: This product is for in vitro research use only