

FOXP1 Polyclonal Antibody

Catalog No: #41621

Package Size: #41621-1 50ul #41621-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FOXP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FOXP1 Polyclonal Antibody detects endogenous levels of FOXP1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FOXP1.
Target Name	FOXP1
Other Names	FOXP1; HSPC215; Forkhead box protein P1
Accession No.	Swiss-Prot: Q9H334NCBI Gene ID: 27086
Uniprot	Q9H334
GeneID	27086;
SDS-PAGE MW	75kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of K562 cells using FOXP1 Polyclonal Antibody

Note: This product is for in vitro research use only