

Neuregulin-2 Polyclonal Antibody

Catalog No: #41623

Package Size: #41623-1 50ul #41623-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Neuregulin-2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Neuregulin-2 Polyclonal Antibody detects endogenous levels of Neuregulin-2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Neuregulin-2.
Target Name	Neuregulin-2
Other Names	NRG2; NTA; Pro-neuregulin-2, membrane-bound isoform; Pro-NRG2
Accession No.	Swiss-Prot: O14511NCBI Gene ID: 9542
Uniprot	O14511
GeneID	9542;
SDS-PAGE MW	90kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HT29 cells using Neuregulin-2 Polyclonal Antibody

Note: This product is for in vitro research use only