

JXCN- cleavage culture

Catalog No: #MSE-3016



Package Size: #MSE-3016-1 10mL/bottle #MSE-3016-2 30mL/bottle

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------|--|
| Product Name | JXCN- cleavage culture |
| Storage | The kit is stored at 2~8°C in the dark and can be stored for 3 months. It is valid for 1 week after opening.Avoid repeated freezing and thawing. |

Product Description

Composition: Purpose: This product is used for invitro culture of mouse fertilized eggs or 2-cell embryos to 8-cell stage. It is a special customized culture medium to support cell division, growth and development of early embryos.Quality Index: 1. pH: 7.2 0.4 (After balance)
2. Osmotic pressure: 30020 mOsm/kg
3. Mouse embryo testMEA:Blastocyst formation rate80% (1-cell was cultured for 96 hours.)
4. Sterilization: filtration method (SAL 10-6)Sample Requirements: Using Protocol: This product must be balanced in a 5% CO2 incubator for at least 2 hours before transferring fertilized eggs or embryos to this culture medium. Embryos will be maintained in this culture medium until the third day. If embryo transferring is scheduled for the third day, it is recommended to transfer the embryo to the uterus with a small amount of this culture medium.Attentions: This product is for research use only.

Note: This product is for in vitro research use only