

mHTF culture medium

Catalog No: #MSE-3018



Package Size: #MSE-3018-1 10mL/bottle #MSE-3018-2 30mL/bottle

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	mHTF culture medium
Storage	The kit is stored at 2~8°C in the dark and can be stored for 3 months. It is valid for 1 week after opening.Avoid repeated freezing and thawing.

Product Description

Composition: Purpose: This product uses HEPES buffer system to maintain the stability of environment (PH value), and is for adjustable medium of all gametes and pre-implantation embryo operation.Quality Index: 1. pH: 7.2 0.4
2. Osmotic pressure: 30020 mOsm/kg
3. Mouse embryo testMEA:Blastocyst formation rate80% (1-cell was cultured for 96 hours.)
4. Sterilization: filtration method (SAL 10-6)Sample Requirements: Using Protocol: The culture medium can only be used in the air environment. There is no need to balance before use.Attentions: This product is for research use only.

Note: This product is for in vitro research use only