

HEXA Polyclonal Antibody

Catalog No: #41641

Package Size: #41641-1 50ul #41641-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HEXA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	HEXA Polyclonal Antibody detects endogenous levels of HEXA protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human HEXA.
Target Name	HEXA
Other Names	HEXA; Beta-hexosaminidase subunit alpha; Beta-N-acetylhexosaminidase subunit alpha; Hexosaminidase subunit A; N-acetyl-beta-glucosaminidase subunit alpha
Accession No.	Swiss-Prot: P06865NCBI Gene ID: 3073
Uniprot	P06865
GeneID	3073;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using HEXA Polyclonal Antibody

Note: This product is for in vitro research use only