

CD107a Polyclonal Antibody

Catalog No: #41667

Package Size: #41667-1 50ul #41667-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD107a Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Rt
Specificity	CD107a Polyclonal Antibody detects endogenous levels of CD107a protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD107a.
Target Name	CD107a
Other Names	LAMP1; Lysosome-associated membrane glycoprotein 1; LAMP-1; Lysosome-associated membrane protein 1; CD107 antigen-like family member A; CD107a
Accession No.	Swiss-Prot: P11279NCBI Gene ID: 3916
Uniprot	P11279
GeneID	3916;
SDS-PAGE MW	120kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of L929 cells using CD107a Polyclonal Antibody

Note: This product is for in vitro research use only