

Nectin 3 Polyclonal Antibody

Catalog No: #41668

Package Size: #41668-1 50ul #41668-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Nectin 3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Nectin 3 Polyclonal Antibody detects endogenous levels of Nectin 3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Nectin 3.
Target Name	Nectin 3
Other Names	PVRL3; PRR3; Poliovirus receptor-related protein 3; CDw113; Nectin-3; CD113
Accession No.	Swiss-Prot: Q9NQS3NCBI Gene ID: 25945
Uniprot	Q9NQS3
GeneID	25945;
SDS-PAGE MW	60kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HeLa cells using Nectin 3 Polyclonal Antibody

Note: This product is for in vitro research use only