

GGT1 Polyclonal Antibody

Catalog No: #41673

Package Size: #41673-1 50ul #41673-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GGT1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	GGT1 Polyclonal Antibody detects endogenous levels of GGT1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human GGT1.
Target Name	GGT1
Other Names	GGT1; GGT; Gamma-glutamyltranspeptidase 1; GGT 1; Gamma-glutamyltransferase 1; Glutathione hydrolase 1; Leukotriene-C4 hydrolase; CD224
Accession No.	Swiss-Prot: P19440NCBI Gene ID: 2678
Uniprot	P19440
GeneID	2678;
SDS-PAGE MW	63kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using GGT1 Polyclonal Antibody

Note: This product is for in vitro research use only