

CD248 Polyclonal Antibody

Catalog No: #41674

Package Size: #41674-1 50ul #41674-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD248 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD248 Polyclonal Antibody detects endogenous levels of CD248 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD248.
Target Name	CD248
Other Names	CD248; CD164L1; TEM1; Endosialin; Tumor endothelial marker 1; CD248
Accession No.	Swiss-Prot: Q9HCU0NCBI Gene ID: 57124
Uniprot	Q9HCU0
GeneID	57124;
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using CD248 Polyclonal Antibody

Note: This product is for in vitro research use only