

CDCP1 Polyclonal Antibody

Catalog No: #41697



Package Size: #41697-1 50ul #41697-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CDCP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CDCP1 Polyclonal Antibody detects endogenous levels of CDCP1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CDCP1.
Target Name	CDCP1
Other Names	CDCP1; TRASK; CUB domain-containing protein 1; Membrane glycoprotein gp140; Subtractive immunization M plus HEP3-associated 135 kDa protein; SIMA135; Transmembrane and associated with src kinases; CD318
Accession No.	Swiss-Prot: Q9H5V8NCBI Gene ID: 64866
Uniprot	Q9H5V8
GeneID	64866;
SDS-PAGE MW	95kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of 293 cells using CDCP1 Polyclonal Antibody

Note: This product is for in vitro research use only