

E2F-4 Polyclonal Antibody

Catalog No: #41716

Package Size: #41716-1 50ul #41716-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	E2F-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	E2F-4 Polyclonal Antibody detects endogenous levels of E2F-4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human E2F-4 around the non-acetylation site of K96.
Target Name	E2F-4
Other Names	E2F4; Transcription factor E2F4; E2F-4
Accession No.	Swiss-Prot: Q16254NCBI Gene ID: 1874
Uniprot	Q16254
GeneID	1874;
SDS-PAGE MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using E2F-4 Polyclonal Antibody

Note: This product is for in vitro research use only