

HSP 90 Polyclonal Antibody

Catalog No: #41721



Package Size: #41721-1 50ul #41721-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HSP 90 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	HSP 90 Polyclonal Antibody detects endogenous levels of HSP 90 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HSP 90 around the non-acetylation site of K292/K284.
Target Name	HSP 90
Other Names	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38; HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
Accession No.	Swiss-Prot: P07900NCBI Gene ID: 3326
Uniprot	P07900
GeneID	3320;
SDS-PAGE MW	85kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using HSP 90 Polyclonal Antibody

Note: This product is for in vitro research use only