

Amphiphysin II Polyclonal Antibody

Catalog No: #41743

Package Size: #41743-1 50ul #41743-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Amphiphysin II Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Amphiphysin II Polyclonal Antibody detects endogenous levels of Amphiphysin II protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Amphiphysin II.
Target Name	Amphiphysin II
Other Names	BIN1; AMPHL; Myc box-dependent-interacting protein 1; Amphiphysin II; Amphiphysin-like protein; Box-dependent myc-interacting protein 1; Bridging integrator 1
Accession No.	Swiss-Prot: O00499NCBI Gene ID: 274
Uniprot	O00499
GeneID	274;
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of 4T1 cells using Amphiphysin II Polyclonal Antibody.

Note: This product is for in vitro research use only