

NTE Polyclonal Antibody

Catalog No: #41751

Package Size: #41751-1 50ul #41751-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NTE Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF(paraffin section)
Species Reactivity	Hu Ms Rt
Specificity	NTE Polyclonal Antibody detects endogenous levels of NTE protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human NTE.
Target Name	NTE
Other Names	PNPLA6; NTE; Neuropathy target esterase; Patatin-like phospholipase domain-containing protein 6
Accession No.	Swiss-Prot: Q8IY17NCBI Gene ID: 10908
Uniprot	Q8IY17
GeneID	10908;
SDS-PAGE MW	150kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


Images




Western Blot analysis of AD293 cells using NTE Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded rat-testis, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded rat-testis, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded rat-testis, antibody was diluted at 1:100

Note: This product is for in vitro research use only