

CRBP II Polyclonal Antibody

Catalog No: #41773

Package Size: #41773-1 50ul #41773-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CRBP II Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CRBP II Polyclonal Antibody detects endogenous levels of CRBP II protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CRBP II.
Target Name	CRBP II
Other Names	RBP2; CRBP2; Retinol-binding protein 2; Cellular retinol-binding protein II; CRBP-II
Accession No.	Swiss-Prot: P50120NCBI Gene ID: 5948
Uniprot	P50120
GeneID	5948;
SDS-PAGE MW	15kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of K562 cells using CRBP II Polyclonal Antibody.

Note: This product is for in vitro research use only