

# Human Sodium- and chloride-dependent taurine transporter (SLC6A6/TAUT) ELISA Kit

Catalog No: #EK7257

Package Size: #EK7257-1 48T #EK7257-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human Sodium- and chloride-dependent taurine transporter (SLC6A6/TAUT) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human (Homo sapiens)
Other Names	MGC10619; MGC131729; TAUT;
Accession No.	P31641
Uniprot	P31641
GeneID	6533;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

## Application Details

Detect Range:0.156-10 ng/mL

Sensitivity:0.056 ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

**Detection Method:**SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate SLC6A6 in samples. An antibody specific for SLC6A6 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anySLC6A6 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SLC6A6 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SLC6A6 bound in the initial step. The color development is stopped and the intensity of the color is measured.

**Product Overview:**Taurine (2-aminoethanesulfonic acid) is a major intracellular amino acid in mammals. It is involved in a number of important physiologic processes, including bile acid conjugation in hepatocytes, modulation of calcium flux and neural excitability, osmoregulation, detoxification, and membrane stabilization. The cells of most organisms respond to hypertonicity by the intracellular accumulation of high concentrations of small organic solutes (osmolytes) that, in contrast to high concentrations of electrolytes, do not perturb the function of macromolecules. The renal medulla is normally the only tissue in mammals that undergoes wide shifts in tonicity. Its hypertonicity when the kidney is excreting a concentrated urine is fundamental to water conservation.

---

Note: This product is for in vitro research use only