

IL-8R β Polyclonal Antibody

Catalog No: #41782

Package Size: #41782-1 50ul #41782-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-8R β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IL-8R β Polyclonal Antibody detects endogenous levels of IL-8R β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human IL-8R β .
Target Name	IL-8R β
Other Names	CXCR2; IL8RB; C-X-C chemokine receptor type 2; CXCR-2; CXCR-2; CDw128b; GRO/MGSA receptor; High affinity interleukin-8 receptor B; IL-8R B; IL-8 receptor type 2; CD182
Accession No.	Swiss-Prot: P25025NCBI Gene ID: 3579
Uniprot	P25025
GeneID	3579;
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of HuvEc cells using IL-8R β Polyclonal Antibody.

Note: This product is for in vitro research use only