

CD163b Polyclonal Antibody

Catalog No: #41791

Package Size: #41791-1 50ul #41791-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD163b Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD163b Polyclonal Antibody detects endogenous levels of CD163b protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD163b.
Target Name	CD163b
Other Names	CD163L1; CD163B; M160; Scavenger receptor cysteine-rich type 1 protein M160; CD163 antigen-like 1; CD163b
Accession No.	Swiss-Prot: Q9NR16NCBI Gene ID: 283316
Uniprot	Q9NR16
GeneID	283316;
SDS-PAGE MW	160kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using CD163b Polyclonal Antibody.

Note: This product is for in vitro research use only