

CD158b2/j Polyclonal Antibody

Catalog No: #41792



Package Size: #41792-1 50ul #41792-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD158b2/j Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD158b2/j Polyclonal Antibody detects endogenous levels of CD158b2/j protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD158b2/j.
Target Name	CD158b2/j
Other Names	KIR2DL3; CD158B2; KIRCL23; NKAT2; Killer cell immunoglobulin-like receptor 2DL3; CD158 antigen-like family member B2; KIR-023GB; Killer inhibitory receptor cl 2-3; MHC class I NK cell receptor; NKAT2a; NKAT2b
Accession No.	Swiss-Prot: P43628/P43631NCBI Gene ID: 3804/100132285
Uniprot	P43628
GeneID	3804;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using CD158b2/j Polyclonal Antibody.

Note: This product is for in vitro research use only