

Inhibin β -A Polyclonal Antibody

Catalog No: #41794



Package Size: #41794-1 50ul #41794-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Inhibin β -A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB, IHC, IF, ELISA
Species Reactivity	Hu Ms Rt
Specificity	Inhibin β -A Polyclonal Antibody detects endogenous levels of Inhibin β -A protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Inhibin β -A.
Target Name	Inhibin β -A
Other Names	INHBA; Inhibin beta A chain; Activin beta-A chain; Erythroid differentiation protein; EDF
Accession No.	Swiss-Prot: P08476NCBI Gene ID: 3624
Uniprot	P08476
GeneID	3624;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-1:2000; IHC: 1:100-1:300; ELISA 1:20000; IF 1:50-200

Images



Western Blot analysis of MCF-7, 293 cells using Inhibin β -A Polyclonal Antibody.

Note: This product is for in vitro research use only