

MDA-7 Polyclonal Antibody

Catalog No: #41796

Package Size: #41796-1 50ul #41796-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MDA-7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	MDA-7 Polyclonal Antibody detects endogenous levels of MDA-7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MDA-7.
Target Name	MDA-7
Other Names	IL24; MDA7; ST16; Interleukin-24; IL-24; Melanoma differentiation-associated gene 7 protein; MDA-7; Suppression of tumorigenicity 16 protein
Accession No.	Swiss-Prot: Q13007NCBI Gene ID: 11009
Uniprot	Q13007
GeneID	11009;
SDS-PAGE MW	24kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using MDA-7 Polyclonal Antibody.

Note: This product is for in vitro research use only