

## TPO Polyclonal Antibody

Catalog No: #41798

Package Size: #41798-1 50ul #41798-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TPO Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	TPO Polyclonal Antibody detects endogenous levels of TPO protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TPO.
Target Name	TPO
Other Names	THPO; MGDF; Thrombopoietin; C-mpl ligand; ML; Megakaryocyte colony-stimulating factor; Megakaryocyte growth and development factor; MGDF; Myeloproliferative leukemia virus oncogene ligand
Accession No.	Swiss-Prot: P40225NCBI Gene ID: 7066
Uniprot	P40225
GeneID	7066;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of MCF-7 cells using TPO Polyclonal Antibody.

---

Note: This product is for in vitro research use only