

PD-ECGF Polyclonal Antibody

Catalog No: #41799

Package Size: #41799-1 50ul #41799-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PD-ECGF Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	PD-ECGF Polyclonal Antibody detects endogenous levels of PD-ECGF protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human PD-ECGF.
Target Name	PD-ECGF
Other Names	TYMP; ECGF1; Thymidine phosphorylase; TP; Gliostatin; Platelet-derived endothelial cell growth factor; PD-ECGF; TdRPase
Accession No.	Swiss-Prot: P19971NCBI Gene ID: 1890
Uniprot	P19971
GeneID	1890;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF-7, 293 cells using PD-ECGF Polyclonal Antibody.

Note: This product is for in vitro research use only