

Neuropilin Polyclonal Antibody

Catalog No: #41826

Package Size: #41826-1 50ul #41826-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neuropilin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Neuropilin Polyclonal Antibody detects endogenous levels of Neuropilin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Neuropilin.
Target Name	Neuropilin
Other Names	NRP1; NRP; VEGF165R; Neuropilin-1; Vascular endothelial cell growth factor 165 receptor; CD304
Accession No.	Swiss-Prot: O14786NCBI Gene ID: 8829
Uniprot	O14786
GeneID	8829;
SDS-PAGE MW	100kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Neuropilin Polyclonal Antibody.

Note: This product is for in vitro research use only