

Integrin $\beta 4$ Polyclonal Antibody

Catalog No: #41840

Package Size: #41840-1 50ul #41840-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin $\beta 4$ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Integrin $\beta 4$ Polyclonal Antibody detects endogenous levels of Integrin $\beta 4$ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Integrin $\beta 4$.
Target Name	Integrin $\beta 4$
Other Names	ITGB4; Integrin beta-4; GP150; CD104
Accession No.	Swiss-Prot: P16144NCBI Gene ID: 3691
Uniprot	P16144
GeneID	3691;
SDS-PAGE MW	200kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of PC12 cells using Integrin $\beta 4$ Polyclonal Antibody.

Note: This product is for in vitro research use only