

CD200 Polyclonal Antibody

Catalog No: #41846

Package Size: #41846-1 50ul #41846-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD200 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD200 Polyclonal Antibody detects endogenous levels of CD200 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD200.
Target Name	CD200
Other Names	CD200; MOX1; MOX2; My033; OX-2 membrane glycoprotein; CD200
Accession No.	Swiss-Prot: P41217NCBI Gene ID: 4345
Uniprot	P41217
GeneID	4345;
SDS-PAGE MW	31kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of PC12 cells using CD200 Polyclonal Antibody.

Note: This product is for in vitro research use only