

Flt3-L Polyclonal Antibody

Catalog No: #41886

Package Size: #41886-1 50ul #41886-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Flt3-L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Flt3-L Polyclonal Antibody detects endogenous levels of Flt3-L protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Flt3-L.
Target Name	Flt3-L
Other Names	FLT3LG; Fms-related tyrosine kinase 3 ligand; Flt3 ligand; Flt3L; SL cytokine
Accession No.	Swiss-Prot: P49771NCBI Gene ID: 2323
Uniprot	P49771
GeneID	2323;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Flt3-L Polyclonal Antibody

Note: This product is for in vitro research use only