

ZNF384 Polyclonal Antibody

Catalog No: #41899



Package Size: #41899-1 50ul #41899-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF384 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Rt
Specificity	ZNF384 Polyclonal Antibody detects endogenous levels of ZNF384 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human ZNF384.
Target Name	ZNF384
Other Names	ZNF384; CAGH1; CIZ; NMP4; TNRC1; Zinc finger protein 384; CAG repeat protein 1; CAS-interacting zinc finger protein; Nuclear matrix transcription factor 4; Nuclear matrix protein 4; Trinucleotide repeat-containing gene 1 protein
Accession No.	Swiss-Prot: Q8TF68NCBI Gene ID: 171017
Uniprot	Q8TF68
GeneID	171017;
SDS-PAGE MW	64kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of VEC cells using ZNF384 Polyclonal Antibody

Note: This product is for in vitro research use only